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HOŘČÍK

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GENERAL	Test name	MAGNESIUM	
	Analysis mode	endpoint mon.	
	Sample type	serum	
	Units	mmol/l	
	Reaction type	increasing	
	Decimals	2	
	Replicates	1	
	Name of assoc. constituent	-	
PROCEDURE	Type of reading	Bichromatic	
	Volumes	Sample	4
		Reagent 1	300
		Reagent 2	-
	Filters	Washing	1.2
		Predilution factor	-
		Main	505
		Reference	600
	Times	Reading 1	312 s
		Reading 2	-
Reagent 2		-	
Postdilution factor		2	
CALIBRATION	Type of calibration	Multiple	
	Calibrator replicates	3	
	Blank replicates	3	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	0.650	
	Kinetic blank limit	-	
	Linearity limit	1.64	

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• LITERATURA

- Chauman UPS and Ray Sarkar BC. Use of calmagite for the determination of traces of magnesium in biological materials. Anal Biochem 1969; 32: 70 – 80.
- Gindler EM and Heth RA. Colorimetric determination with bound "calmagite" of magnesium in human blood serum. Clin Chem 1971; 17: 662.
- Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
- Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

BioSystems S. A.
Costa Brava 30
Barcelona, SPAIN

Distibutor v ČR: dot®diagnostics, s.r.o.

Ruzyňská 519/16
CZ - 161 00 Praha 6
Tel.: +420 235 318 612
Fax: +420 235 318 614
e-mail: dotdiag@dotdiag.cz