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## ALP (AMP) IFCC

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<b>GENERAL</b>	Test name	ALP-AMP	
	Analysis mode	kinetic mon.	
	Sample type	serum	
	Units	μkat/l	
	Reaction type	increasing	
	Decimals	2	
	Replicates	1	
Name of assoc. constituent		-	
<b>PROCEDURE</b>	Type of reading	monoch.	
	Volumes	Sample	6
		Reagent 1	300
		Reagent 2	-
		Washing	1.2
	Filters	Predilution factor	-
		Main	405
	Times	Reference	-
		Reading 1	72 s
		Reading 2	216 s
Reagent 2		-	
Postdilution factor		2	
<b>CALIBRATION</b>	Type of calibration	Multiple	
	Calibrator replicates	3	
	Blank replicates	3	
	Calibration curve	-	
<b>OPTIONS</b>	Blank absorbance limit	0.600	
	Kinetic blank limit	-	
	Linearity limit	20	

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## • LITERATURA

- IFCC methods for the measurement of catalytic concentration of enzymes. Part 5: IFCC methods for alkaline phosphatase. *J Clin Chem Clin Biochem* 1983; 21: 731 – 748.
- Rosalki SB, Foo AY, Burlina A, et al. Multicenter evaluation of iso-ALP test kit for measurement of bone alkaline phosphatase activity in serum and plasma. *Clin Chem* 1993; 39: 648 – 652.
- Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.

BioSystems S. A.  
Costa Brava 30  
Barcelona, SPAIN

Distibutor v ČR: dot®diagnostics, s.r.o.  
Ruzyňská 519/16  
CZ - 161 00 Praha 6  
Tel.: +420 235 318 612  
Fax: +420 235 318 614  
e-mail: dotdiag@dotdiag.cz