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KYSELINA MOČOVÁ

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GENERAL	Test name	URIC ACID
	Analysis mode	endpoint mon.
	Sample type	serum
	Units	$\mu\text{mol/l}$
	Reaction type	increasing
	Decimals	0
	Replicates	1
Name of assoc. constituent		-
PROCEDURE	Type of reading	bichrom.
	Volumes	
	Sample	7.5
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	505
	Reference	670
Times	Reading 1	312 s
	Reading 2	-
	Reagent 2	-
	Postdilution factor	2
CALIBRATION	Type of calibration	Multiple
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.200
	Kinetic blank limit	-
	Linearity limit	1487 $\mu\text{mol/l}$

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
• LITERATURA

- Barham D, Trinder P. An improved colour reagent for the determination of blood glucose by oxidase system. *Analyst* 1972; 27: 142 – 145.
- Fossati P, Prencipe L, Berti G. Use of 3,5-dichloro-2-hydroxybenzenesulfonicacid/4-aminophenazone chromogenic system in direct enzymic assay of uric acid in serum and urine. *Clin Chem* 1980; 26: 227 – 231.
- Tietz NW, ed. Clinical guide to laboratory test, 2nd ed. Saunders Co, 1991.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

BioSystems S. A.
Costa Brava 30
Barcelona, SPAIN

Distibutor v ČR: dot®diagnostics, s.r.o.

Ruzyňská 519/16
CZ - 161 00 Praha 6
Tel.: +420 235 318 612
Fax: +420 235 318 614
e-mail: dotdiag@dotdiag.cz