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TRIGLYCERIDY

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GENERAL	Test name	TRIGLYCERIDES
	Analysis mode	endpoint mon.
	Sample type	serum
	Units	mmol/l
	Reaction type	increasing
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	bichrom.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	505
	Reference	670
Times	Reading 1	312 s
	Reading 2	-
	Reagent 2	-
	Postdilution factor	2
CALIBRATION	Type of calibration	multiple
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.150
	Kinetic blank limit	-
	Linearity limit	6.78

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• LITERATURA

- Bucolo G and David H. Quantitative determination of serum triglycerides by use of enzymes. *Clin Chem* 1973; 19: 476 – 482.
- Fossati P and Prencipe L. Serum triglycerides determined colorimetrically with an enzyme that produces hydrogen peroxide. *Clin Chem* 1982; 28: 2077 – 2080.
- National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol education Program (NCEP Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: National Heart, Lung, and Blood Institute; 2001.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

BioSystems S. A.
Costa Brava 30
Barcelona, SPAIN

Distibutor v ČR: dot®diagnostics, s.r.o.

Ruzyňská 519/16
CZ - 161 00 Praha 6
Tel.: +420 235 318 612
Fax: +420 235 318 614
e-mail: dotdiag@dotdiag.cz