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## α-AMYLÁZA

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<b>GENERAL</b>	Test name	AMYLASE
	Analysis mode	kinetic mon.
	Sample type	Serum
	Units	μkat/l
	Reaction type	increasing
	Decimals	2
	Replicates	1
Name of assoc. constituent		-
<b>PROCEDURE</b>	Type of reading	monoch.
	Volumes	
	Sample	6
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	405
	Reference	-
Times	Reading 1	72 s
	Reading 2	216 s
	Reagent 2	-
	Postdilution factor	2
<b>CALIBRATION</b>	Type of calibration	Multiple
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
<b>OPTIONS</b>	Blank absorbance limit	0.200
	Kinetic blank limit	-
	Linearity limit	21.7

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• **LITERATURA**

- Winn-Deen ES, David H, Sigler G and Chavez R. Development of a direct assay for α-amylase. *Clin Chem* 1988; 34: 2005 – 2008.
- Gella FJ, Gubern G, Vidal R, Canalias F. Determination of total and pancreatic α-amylase in human serum with 2-chloro-4-nitrophenyl-α-D-maltotriose as substrate. *Clin Chem Acta* 1997; 259: 147 – 160.
- Gubern G, Balsells D, Ferragut R, Galán A, Gella FJ, et al. Procedimiento recomendado para la determinación en rutina de la concentración catalítica de α-amilasa en suero sanguíneo humano. *Quim Clin* 1996; 15: 51 – 52.
- Balsells D, Gella FJ, Gubern G, Canalias F. Reference values for α-amylase in human serum and urine using 2-chloro-4-nitrophenyl-α-D-maltotriose as substrate. *Clin Chem Acta* 1998; 274: 213 – 217.
- Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CS, Ashood ER. WB Saunders Co., 1999.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

**BioSystems S. A.**  
Costa Brava 30  
Barcelona, SPAIN

**Distibutor v ČR: dot®diagnostics, s.r.o.**

Ruzyňská 519/16  
CZ - 161 00 Praha 6  
Tel.: +420 235 318 612  
Fax: +420 235 318 614  
e-mail: [dotdiag@dotdiag.cz](mailto:dotdiag@dotdiag.cz)