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HDL - CHOLESTEROL

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GENERAL	Test name	CHOLESTEROL HDL DIRECT	
	Analysis mode	bir.differential	
	Sample type	serum	
	Units	mmol/l	
	Reaction type	increasing	
	Decimals	2	
	Replicates	1	
	Name of assoc. constituent	-	
	PROCEDURE	Type of reading	monoch.
		Volumes	
Sample		3	
Reagent 1		300	
Reagent 2		100	
Washing		1.2	
Predilution factor		-	
Filters			
Main		535	
Reference		-	
Times	Reading 1	168 s	
	Reading 2	480 s	
	Reagent 2	192 s	
	Postdilution factor	2	
CALIBRATION	Type of calibration	multiple	
	Calibrator replicates	3	
	Blank replicates	3	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	0.05	
	Kinetic blank limit	-	
	Linearity limit	5.18	

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● LITERATURA

- Warnick GR Nauck M, Rifai N. Evolution of methods for measurement of HDL-cholesterol: from ultracentrifugation to homogeneous assays. *Clin Chem* 2001; 47:1579 - 1596.
- National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol education Program (NCEP Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: National Heart, Lung, and Blood Institute; 2001.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Tietz NW. Clinical guide to laboratory tests, 2nd ed. Saunders Co, 1991.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

BioSystems S. A.
Costa Brava 30
Barcelona, SPAIN

Distibutor v ČR: dot®diagnostics, s.r.o.

Ruzyňská 519/16
CZ - 161 00 Praha 6
Tel.: +420 235 318 612
Fax: +420 235 318 614
e-mail: dotdiag@dotdiag.cz