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GENERAL	Test name	PROTEIN URINE
	Analysis mode	endpoint mon.
	Sample type	urine
	Units	g/l
	Reaction type	increasing
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
PROCEDURE	Type of reading	bichrom.
	Volumes	
	Sample	6
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	600
	Reference	670
Times	Reading 1	312 s
	Reading 2	-
	Reagent 2	-
	Postdilution factor	2
	CALIBRATION	Type of calibration
Calibrator replicates		3
Blank replicates		3
Calibration curve		-
OPTIONS	Blank absorbance limit	0.120
	Kinetic blank limit	-
	Linearity limit	4

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• LITERATURA

- Watanabe N et al. Urinary Protein as measured with a pyrogallol red-molybdate complex, manually and in a Hitachi 726 automated analyzer. *Clin Chem* 1986; 32: 1551 – 1544.
- Orsonneau JL et al. An improved pyrogallol red-molybdate method for the determining total urinary protein. *Clin Chem* 1989; 35: 2233 – 2236.
- Tietz NW et al. Clinical guide to laboratory. 3rd ed Saunders Co, 1999.
- Young DS. Effects of drugs on clinical laboratory tests 4th ed. AACC Press, 1995.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

BioSystems S. A.
Costa Brava 30
Barcelona, SPAIN

Distibutor v ČR: dot®diagnostics, s.r.o.
Ruzyňská 519/16
CZ - 161 00 Praha 6
Tel.: +420 235 318 612
Fax: +420 235 318 614
e-mail: dotdiag@dotdiag.cz