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LDL- CHOLESTEROL

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GENERAL	Test name	CHOLESTEROL LDL DIRECT
	Analysis mode	differential bir.
	Sample type	serum
	Units	mmol/l
	Reaction type	increasing
	Decimals	2
	Replicates	1
Name of assoc. constituent	-	
PROCEDURE	Type of reading	monoch.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	100
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	535
	Reference	-
Times	Reading 1	288 s
	Reading 2	480 s
	Reagent 2	312 s
	Postdilution factor	2
CALIBRATION	Type of calibration	multiple
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.100
	Kinetic blank limit	-
	Linearity limit	25.6

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
● **LITERATURA**

- Nauck M, Warnick GR, Rifai N. Methods for measurement of LDL-cholesterol: a critical assessment of direct measurement by homogeneous assays versus calculation. *Clin Chem* 2002; 48: 236 – 254.
- National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP Expert Panel on Detection, Evaluation, and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: National Heart, Lung, and Blood Institute; 2001.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Tietz NW. Clinical guide to laboratory tests, 3rd ed. Saunders Co., 1999.
- Friedman and Young. Effects of disease on clinical laboratory test, 3th ed. AACC Press, 1997.

Distibutor v ČR: dot®diagnostics, s.r.o.

Ruzyňská 519/16
CZ - 161 00 Praha 6
Tel.: +420 235 318 612
Fax: +420 235 318 614
e-mail: dotdiag@dotdiag.cz