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Imunoglobulin G

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Jako vzorek použijte **sérum** nebo **plasmu** (EDTA, heparin).

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GENERAL	Test name	IgG
	Analysis mode	endpoint mon.
	Sample type	serum
	Units	g/l
	Reaction type	increasing
	Turbidim. test	yes
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
	PROCEDURE	Type of reading
Volumes		
Sample		3
Reagent 1		440
Reagent 2		-
Washing		1.II
Predilution factor		-
Filters		
Main		535
Reference		-
Times	Reading 1	312 s
	Reading 2	-
	Reagent 2	-
	Postdilution factor	1.IV
CALIBRATION	Type of calibration	multiple
	Calibrator number	5
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	increasing polygonal
OPTIONS	Blank absorbance limit	0.300
	Kinetic blank limit	-
	Linearity limit	-

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
• LITERATURA

- Narayanan S. Method-comparison studies on immunoglobulins. *Clin Chem* 1982; 28: 1528 – 1531.
- Price CP, Spencer K and Whicher J. Light-scattering immunoassay of specific proteins: a review. *Ann Clin Biochem* 1983; 20: 1 – 14.
- Dati F et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference range for 14 proteins in serum based on the standardization against the IFCC/CAP reference material (CRM470). *Eur J Clin Chem Clin Biochem* 1996; 34: 517 – 520.
- Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th, ed. AACC Press, 1997.

BioSystems S. A.
Costa Brava 30
Barcelona, SPAIN

Distibutor v ČR: dot®diagnostics, s.r.o.
Ruzyňská 519/16
CZ - 161 00 Praha 6
Tel.: +420 235 318 612
Fax: +420 235 318 614
e-mail: dotdiag@dotdiag.cz