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## Imunoglobulin M

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Jako vzorek použijte **sérum** nebo **plasmu** (EDTA, heparin).

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GENERAL	Test name	IgM
	Analysis mode	endpoint mon.
	Sample type	serum
	Units	g/l
	Reaction type	increasing
	Turbidim. test	yes
	Decimals	2
	Replicates	1
	Name of assoc. constituent	-
	PROCEDURE	Type of reading
Volumes		
Sample		5
Reagent 1		300
Reagent 2		-
Washing		1.II
Predilution factor		-
Filters		
Main		340
Reference		670
Times	Reading 1	600 s
	Reading 2	-
	Reagent 2	-
	Postdilution factor	2
	CALIBRATION	Type of calibration
Calibrator number		5
Calibrator replicates		3
Blank replicates		3
Calibration curve		increasing polygonal
OPTIONS	Blank absorbance limit	0.300
	Kinetic blank limit	-
	Linearity limit	-
	Linearity limit	-

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● **LITERATURA**

- Narayanan S. Method-comparison studies on immunoglobulins. *Clin Chem* 1982; 28: 1528 – 1531.
- Price CP, Spencer K and Whicher J. Light-scattering immunoassay of specific proteins: a review. *Ann Clin Biochem* 1983; 20: 1 – 14.
- Dati F et al. Consensus of a group of professional societies and diagnostic companies on guidelines for interim reference range for 14 proteins in serum based on the standardization against the IFCC/CAP reference material (CRM470). *Eur J Clin Chem Clin Biochem* 1996; 34: 517 – 520.
- Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th, ed. AACC Press, 1997.

**Distibutor v ČR: dot®diagnostics, s.r.o.**  
Ruzyňská 519/16  
CZ - 161 00 Praha 6  
Tel.: +420 235 318 612  
Fax: +420 235 318 614  
e-mail: [dotdiag@dotdiag.cz](mailto:dotdiag@dotdiag.cz)