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CRP

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GENERAL	Test name	CRP	
	Analysis mode	endpoint mon.	
	Sample type	serum	
	Units	mg/l	
	Reaction type	increasing	
	Turbidim. test	yes	
	Decimals	1	
	Replicates	1	
	Name of assoc. constituent	-	
	PROCEDURE	Type of reading	bichrom.
Volumes		Sample	3
		Reagent 1	440
		Reagent 2	-
		Washing	1.2
		Predilution factor	-
Filters		Main	535
		Reference	670
Times		Reading 1	192 s
		Reading 2	-
	Reagent 2	-	
	Postdilution factor	1.4	
CALIBRATION	Type of calibration	Specific	
	Calibrator number	1	
	Calibrator replicates	3	
	Blank replicates	3	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	0.700	
	Kinetic blank limit	-	
	Linearity limit	150	

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• LITERATURA

- Kindmark C-O. The concentration of C-Reactive Protein in sera from healthy individuals. *Scand J Clin Lab Invest* 1972; 29: 407 – 411.
- Grange J, Roch AM, Quasch GA. Nephelometric assay of antigens and antibodies with latex particles. *J Immunol Methods* 1977; 18: 365 – 375.
- Price CP, Trull AK, Berry D, Gorman EG. Development and validation of a particle-enhanced turbidimetric immunoassay for C-reactive protein. *J Immunol Methods* 1987; 99: 205 – 211.
- Otsuji S, Shibata H, Umeda M. Turbidimetric immunoassay of serum C-reactive protein. *Clin Chem* 1982; 28: 2121 – 4.
- Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1999.
- Young DS. Effects of drugs on clinical laboratory, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

BioSystems S. A.
Costa Brava 30
Barcelona, SPAIN

Distibutor v ČR: dot®diagnostics, s.r.o.

Ruzyňská 519/16
CZ - 161 00 Praha 6
Tel.: +420 235 318 612
Fax: +420 235 318 614
e-mail: dotdiag@dotdiag.cz