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## ANTISTREPTOLYSIN O

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|                    |                            |                    |          |
|--------------------|----------------------------|--------------------|----------|
| <b>GENERAL</b>     | Test name                  | ASO                |          |
|                    | Analysis mode              | endpoint mon.      |          |
|                    | Sample type                | serum              |          |
|                    | Units                      | IU/ml              |          |
|                    | Reaction type              | increasing         |          |
|                    | Turbidim. test             | yes                |          |
|                    | Decimals                   | 0                  |          |
|                    | Replicates                 | 1                  |          |
|                    | Name of assoc. constituent | -                  |          |
|                    | <b>PROCEDURE</b>           | Type of reading    | bichrom. |
| Volumes            |                            | Sample             | 3        |
|                    |                            | Reagent 1          | 300      |
|                    |                            | Reagent 2          | -        |
|                    |                            | Washing            | 1.2      |
|                    |                            | Predilution factor | -        |
| Filters            |                            | Main               | 535      |
|                    |                            | Reference          | 670      |
| Times              |                            | Reading 1          | 192 s    |
|                    |                            | Reading 2          | -        |
|                    | Reagent 2                  | -                  |          |
|                    | Postdilution factor        | 2                  |          |
| <b>CALIBRATION</b> | Type of calibration        | specific           |          |
|                    | Calibrator number          | 1                  |          |
|                    | Calibrator replicates      | 3                  |          |
|                    | Blank replicates           | 3                  |          |
|                    | Calibration curve          | -                  |          |
| <b>OPTIONS</b>     | Blank absorbance limit     | 0.700              |          |
|                    | Kinetic blank limit        | -                  |          |
|                    | Linearity limit            | 800                |          |

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## • LITERATURA

- Borque L, Rus A, Dubois H. Automated determination of streptolysin O antibodies by turbidimetric latex immunoassay method. *J Clin Immunoassay* 1992; 15:182 - 6.
- Klein GC, Baker CN, Jones WL. Upper limits of normal antistreptolysin O and antideoxyribonuclease B titers. *Appl Microbiol* 1971; 21: 758 - 60.
- Bisno AL. Group A streptococcal infections and acute rheumatic fever. *N Engl J Med* 1991; 325: 783 - 93.
- Stevens DL. Invasive group A streptococcal disease. *Clin Infect Dis* 1992; 14: 2 -11.
- Immunology and Serology in Laboratory Medicine, 2<sup>nd</sup> edition. Turgeon ML. Mosby, 1996.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACCC Press, 1997.
- Young DS. Effects of drugs on clinical laboratory, 3th ed. AACCC Press, 1997.

**Distibutor v ČR: dot®diagnostics, s.r.o.**

Ruzyňská 519/16  
CZ - 161 00 Praha 6  
Tel.: +420 235 318 612  
Fax: +420 235 318 614  
e-mail: dotdiag@dotdiag.cz