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GLUKÓZA GOD

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GENERAL	Test name	GLUCOSE	
	Analysis mode	endpoint mon.	
	Sample type	serum	
	Units	mmol/l	
	Reaction type	increasing	
	Decimals	2	
	Replicates	1	
	Name of assoc. constituent	-	
PROCEDURE	Type of reading	bichromatic	
	Volumes	Sample	3
		Reagent 1	300
		Reagent 2	-
	Filters	Washing	1.2
		Predilution factor	-
		Main	505
	Times	Reference	670
		Reading 1	312 s
		Reading 2	-
Reagent 2		-	
	Postdilution factor	2	
CALIBRATION	Type of calibration	Multiple	
	Calibrator replicates	3	
	Blank replicates	3	
	Calibration curve	-	
OPTIONS	Blank absorbance limit	0.150	
	Kinetic blank limit	-	
	Linearity limit	27.5	

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• LITERATURA

- Trinder P. Determination of glucose in blood using glucose oxidase with an alternative oxygen acceptor. Ann Clin Biochem 1966; 6: 24-27
- Tietz NW. Clinical guide to laboratory tests, 2nd ed. Saunders Co, 1991.
- National Diabetes Data Group: Classification and diagnosis of diabetes mellitus and other categories of glucose intolerance. Diabetes 1979; 28: 1039-1057.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3rd ed. AACC Press, 1997.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.

BioSystems S. A.
Costa Brava 30
Barcelona, SPAIN

Distibutor v ČR: dot®diagnostics, s.r.o.
Ruzyňská 519/16
CZ - 161 00 Praha 6
Tel.: +420 235 318 612
Fax: +420 235 318 614
e-mail: dotdiag@dotdiag.cz