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## CHOLESTEROL

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GENERAL	Test name	CHOLESTEROL
	Analysis mode	endpoint mon.
	Sample type	serum
	Units	mmol/l
	Reaction type	increasing
	Decimals	2
	Replicates	1
Name of assoc. constituent		-
PROCEDURE	Type of reading	bichrom.
	Volumes	
	Sample	3
	Reagent 1	300
	Reagent 2	-
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	505
	Reference	670
	Times	
	Reading 1	312 s
Reading 2	-	
Reagent 2	-	
Postdilution factor	2	
CALIBRATION	Type of calibration	Multiple
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.200
	Kinetic blank limit	-
	Linearity limit	26

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
## • LITERATURA

- Allain CC Poon LS, Chan CSG, Richmond W and Fu PC. Enzymatic determination of total serum cholesterol. *Clin. Chem* 1974; 20: 470 - 475.
- Meiatini F, Prencipe L, Bardelli F, Giannini G and Tarli P. the 4-hydroxybenzoate/4-aminophenazone chromogenic system used in the enzymic determination of serum cholesterol. *Clin Chem* 1978; 24: 2161 - 2165.
- National Cholesterol Education Program Expert Panel. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (ATP III). NIH Publication. Bethesda: National Heart, Lung, and Blood Institute; 2001.
- Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACC Press, 1995.
- Tietz NW. Clinical guide to laboratory tests, 2nd ed. Saunders Co, 1991.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

**BioSystems S. A.**  
Costa Brava 30  
Barcelona, SPAIN

**Distibutor v ČR: dot®diagnostics, s.r.o.**

Ruzyňská 519/16  
CZ - 161 00 Praha 6  
Tel.: +420 235 318 612  
Fax: +420 235 318 614  
e-mail: [dotdiag@dotdiag.cz](mailto:dotdiag@dotdiag.cz)