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FOSFOR

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Moč sbírejte do 10ml 10% HCl.

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GENERAL	Test name	PHOSPHORUS
	Analysis mode	differential bir.
	Sample type	serum
	Units	mmol/l
	Reaction type	increasing
	Decimals	2
	Replicates	1
Name of assoc. constituent		-
PROCEDURE	Type of reading	monoch.
	Volumes	
	Sample	3
	Reagent 1	210
	Reagent 2	90
	Washing	1.2
	Predilution factor	-
	Filters	
	Main	340
	Reference	-
Times	Reading 1	72 s
	Reading 2	312 s
	Reagent 2	96 s
	Postdilution factor	2
CALIBRATION	Type of calibration	Multiple
	Calibrator replicates	3
	Blank replicates	3
	Calibration curve	-
OPTIONS	Blank absorbance limit	0.500
	Kinetic blank limit	-
	Linearity limit	6.46

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• LITERATURA

- Gamst O and Try K. Determination of serum-phosphate without deproteinization by ultraviolet spectrophotometry of the phosphomolybdc acid complex. *Scand J Clin Lab Invest* 1980; 40: 483 – 486.
- Muñoz MA, Balón M and Fernández C. direct determination of inorganic phosphorus in serum with a single reagent. *Clin Chem* 1983; 29: 372 – 374.
- Tietz Textbook of Clinical Chemistry, 3rd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1999.
- Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

BioSystems S. A.
Costa Brava 30
Barcelona, SPAIN

Distibutor v ČR: dot[®]diagnostics, s.r.o.
Ruzyňská 519/16
CZ - 161 00 Praha 6
Tel.: +420 235 318 612
Fax: +420 235 318 614
e-mail: dotdiag@dotdiag.cz