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## ŽELEZO FERROZIN

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<b>GENERAL</b>	Test name	IRON FERROZINE	
	Analysis mode	differential bir.	
	Sample type	serum	
	Units	µmol/l	
	Reaction type	increasing	
	Decimals	1	
	Replicates	1	
	Name of assoc. constituent	-	
	<b>PROCEDURE</b>	Type of reading	monoch.
		Volumes	
Sample		40	
Reagent 1		240	
Reagent 2		60	
Washing		1.2	
Predilution factor		-	
Filters			
Main		560	
Reference		-	
<b>Times</b>	Reading 1	72 s	
	Reading 2	408 s	
	Reagent 2	96 s	
	Postdilution factor	1,4	
	<b>CALIBRATION</b>	Type of calibration	Multiple
Calibrator replicates		3	
Blank replicates		3	
Calibration curve		-	
<b>OPTIONS</b>	Blank absorbance limit	0.050	
	Kinetic blank limit	-	
	Linearity limit	179	

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• **LITERATURA**

- Stookey LL. Ferrozine-A next spectrophotometric reagent for iron. Anal Chem 1970; 42: 779 – 81.
- Itano M. Serum Iron Survey. Am J Clin Pathol 1978; 70: 516 – 522.
- Artiss JD, Vinogradov S, Zak B. Spectrophotometric study of several sensitive reagents for serum iron. Clin Biochem 1981; 14: 311 – 315.
- Tietz Textbook of Clinical Chemistry, 2nd edition. Burtis CA, Ashwood ER. WB Saunders Co., 1994.
- Young DS. Effects of drugs on clinical laboratory tests, 3th ed. AACC Press, 1997.
- Friedman and Young. Effects of disease on clinical laboratory tests, 3th ed. AACC Press, 1997.

**BioSystems S. A.**  
Costa Brava 30  
Barcelona, SPAIN

**Distibutor v ČR: dot®diagnostics, s.r.o.**

Ruzyňská 519/16  
CZ - 161 00 Praha 6  
Tel.: +420 235 318 612  
Fax: +420 235 318 614  
e-mail: [dotdiag@dotdiag.cz](mailto:dotdiag@dotdiag.cz)